

## Product datasheet for **RC216820L3V**

### OS9 (NM\_001017958) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	OS9 (NM_001017958) Human Tagged ORF Clone Lentiviral Particle
Symbol:	OS9
Synonyms:	ERLEC2; OS-9
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001017958
ORF Size:	1956 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC216820).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001017958.1</a> , <a href="#">NP_001017958.1</a>
RefSeq Size:	2676 bp
RefSeq ORF:	1959 bp
Locus ID:	10956
UniProt ID:	<a href="#">Q13438</a>
Cytogenetics:	12q13.3-q14.1
Protein Families:	Transmembrane
MW:	73.82 kDa



[View online »](#)

**Gene Summary:**

This gene encodes a protein that is highly expressed in osteosarcomas. This protein binds to the hypoxia-inducible factor 1 (HIF-1), a key regulator of the hypoxic response and angiogenesis, and promotes the degradation of one of its subunits. Alternate transcriptional splice variants, encoding different isoforms, have been characterized. [provided by RefSeq, Jul 2008]