

Product datasheet for RC216830

OR4A16 (NM_001005274) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	OR4A16 (NM_001005274) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	OR4A16
Synonyms:	OR4A16Q; OR11-117
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>RC216830 representing NM_001005274 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGAGACCAAGCAGCAATGTTACAGAATTTGTCCTCTGGGCCTCACTCAAGATCCTGATGTGAAAAAA
CATTATTTGTCATGTTTTACTCATATACATTGTGACAATGGTGGGAAACCTCCTCATTGGGTGACTAC
TATTGGCAGCCCCTCCTGGGCTCCCTAATGTACTTCTCCTTGCCACTTGTCACTTATGGATGCCATA
TATCCACTGCCATGTCACCCAAATTGATGATAGACTTACTCTGTGATAAAATCGCTATTTCTTGTGAG
CTTGATGGGTGAGTCTTCATAGAACACTTACTTGGTGGTGCAGAGGTCTTCTTTGGTGGTGTGGC
CTATGATCGCTATGTGGCTATCTCTAAGCCGCTGCACTATTTGAACATCATGAATCGACTGGTTGCATC
CTTCTGTTGGTGGTGGCCATGATTGGAGGTTTTGTGCACTCTGTGGTTCAAATTGTCTTTCTGTACAGTC
TACCAATCTGTGGCCCAATGTTATTGACCACTCTGTCTGTGACATGTACCCATTGTTGGAAGTGTGTG
CCTTGACACCTACTTTATAGGACTCACTGTGGTTGCCAATGGTGGAAATAATTTGTATGGTCATCTTACC
TTTCTGCTAATCTCCTGTGGAGTCATCCTAAACTTCTTAAACTTACAGTCAGGAAGAGAGGCATAAAG
CCCTGCCTACCTGCATCTCCACATCATTGTGGTTGCCCTCGTTTTTGTCCCTGATTTTTATGTATGT
TAGACCCGTTTCCAACCTTCCCTTTGATAAATTAATGACTGTGTTTTATTCAATTATCACACTCATGTTG
AATCCTTTAATACTCGTTGAGACAATCAGAGATGAAAAATGCTATGAAAAATCTGGTGTGAAAAAGT
TAAGTATAGTTAGAAAAAGAGTATCTCCACACTGAACATATTTATTCTAGTTCTAAGGCAACAAATAG
GCGG

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >RC216830 representing NM_001005274
Red=Cloning site Green=Tags(s)

MRPSSNVTEFVLLGLTQDPDVKKTLFVMFLLIYIVTMVGNLLIWVTTIGSPSLGSLMYFFLAYLSLMDAI
 YSTAMSPKLMIDLCDKIAIISLSACMGQLFIEHLLGGAEVFLLVVMAYDRYVAISKPLHYLNIMNRLVCI
 LLLVVAMIGGFVHSVVQIVFLYSLPICGPNVIDHSVCDMYLLELLCLDTYFIFGLTVVANGGIICMVIFT
 FLLISCGVILNFKTYSQEERHKALPTCISHIIVVALVFVPCIFMYVYRVPVSNFPFDKLMTVFYSIITLML
 NPLIYSLRQSEMKNAMKNLWCEKLSIVRKRVSPTLNIFIPSSKATNRR

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mk8003_f03.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001005274

ORF Size: 984 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001005274.1](#), [NP_001005274.1](#)

RefSeq Size: 987 bp

RefSeq ORF: 987 bp

Locus ID: 81327

UniProt ID: [Q8NH70](#)

Cytogenetics: 11q11


Protein Families: Transmembrane

Protein Pathways: Olfactory transduction

MW: 36.7 kDa

Gene Summary: Olfactory receptors interact with odorant molecules in the nose, to initiate a neuronal response that triggers the perception of a smell. The olfactory receptor proteins are members of a large family of G-protein-coupled receptors (GPCR) arising from single coding-exon genes. Olfactory receptors share a 7-transmembrane domain structure with many neurotransmitter and hormone receptors and are responsible for the recognition and G protein-mediated transduction of odorant signals. The olfactory receptor gene family is the largest in the genome. The nomenclature assigned to the olfactory receptor genes and proteins for this organism is independent of other organisms. [provided by RefSeq, Jul 2008]

Product images:



Circular map for RC216830