

Product datasheet for **RC216837L3V**

zinc finger protein 655 (ZNF655) (NM_001083956) Human Tagged ORF Clone Lentiviral Particle

Product data:

| | |
|---------------------------|--|
| Product Type: | Lentiviral Particles |
| Product Name: | zinc finger protein 655 (ZNF655) (NM_001083956) Human Tagged ORF Clone Lentiviral Particle |
| Symbol: | zinc finger protein 655 |
| Synonyms: | VIK; VIK-1 |
| Mammalian Cell Selection: | Puromycin |
| Vector: | pLenti-C-Myc-DDK-P2A-Puro (PS100092) |
| Tag: | Myc-DDK |
| ACCN: | NM_001083956 |
| ORF Size: | 1578 bp |
| ORF Nucleotide Sequence: | The ORF insert of this clone is exactly the same as(RC216837). |
| OTI Disclaimer: | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info |
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| RefSeq: | NM_001083956.1 , NP_001077425.1 |
| RefSeq Size: | 4863 bp |
| RefSeq ORF: | 1581 bp |
| Locus ID: | 79027 |
| UniProt ID: | Q8N720 |
| Cytogenetics: | 7q22.1 |
| Protein Families: | Transcription Factors |



[View online »](#)

MW: 61.4 kDa

Gene Summary: This gene encodes a zinc finger protein. The zinc finger proteins are involved in DNA binding and protein-protein interactions. Multiple alternatively spliced transcript variants encoding distinct isoforms have been found for this gene. [provided by RefSeq, Jul 2008]