

## Product datasheet for **RC216856L3V**

### **NIPP1 (PPP1R8) (NM\_002713) Human Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	NIPP1 (PPP1R8) (NM_002713) Human Tagged ORF Clone Lentiviral Particle
Symbol:	NIPP1
Synonyms:	ARD-1; ARD1; NIPP-1; NIPP1; PRO2047
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_002713
ORF Size:	381 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC216856).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_002713.3</a>
RefSeq Size:	2558 bp
RefSeq ORF:	384 bp
Locus ID:	5511
UniProt ID:	<a href="#">Q12972</a>
Cytogenetics:	1p35.3
Protein Families:	Druggable Genome, Transcription Factors
MW:	13.3 kDa


[View online »](#)

**Gene Summary:**

This gene, through alternative splicing, encodes three different isoforms. Two of the protein isoforms encoded by this gene are specific inhibitors of type 1 serine/threonine protein phosphatases and can bind but not cleave RNA. The third protein isoform lacks the phosphatase inhibitory function but is a single-strand endoribonuclease comparable to RNase E of *E. coli*. This isoform requires magnesium for its function and cleaves specific sites in A+U-rich regions of RNA. [provided by RefSeq, Jul 2008]