

Product datasheet for **RC216896**

SDF4 (NM_016547) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: SDF4 (NM_016547) Human Tagged ORF Clone
Tag: Myc-DDK
Symbol: SDF4
Synonyms: Cab45; SDF-4
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RC216896 representing NM_016547
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGC**C

ATGGTCTGGCCCTGGGTGGCGATGGCGTCCAGGTGGGGTCCCCTCATTGGCCTGGCTCCGTGCTGCCTCT
GGCTCCTGGGGCAGTCCTTCTGATGGACGCGTCTGCACGGCCTGCCAACCCTCGTCCACTCGAGAGAG
AGTAGCCAACAGGGAGGAGAATGAGATCCTGCCCCAGACCACCTGAACGGGGTGAAGCTGGAGATGGAC
GGGCACCTCAATCGCGGCTTCCACCAGGAGGTCTCCTAGGCAAGGACCTGGGTGGCTTTGATGAGGACG
CGGAGCCGCGGGAGCCGGAGGAAGCTGATGGTCATCTTTTCCAAGGTGGATGTGAACACTGACCGGAA
GATCAGTGCCAAGGAGATGCAGCGCTGGATCATGGAGAAGACGGCCGAGCACTTCCAGGAGGCCATGGAG
GAGAGCAAGACACACTTCCGCGCCGTGGACCCTGACGGGGACGGTCACGTGTCTTGGGACGAGTATAAGG
TGAAGTTTTTGGCGAGTAAAGGCCATAGCGAGAAGGAGTTGCCGACGCCATCAGGCTCAACGAGGAACT
CAAAGTGGATGAGGAAACACAGGAAGTCTGGAGAACCTGAAGGACCGCTGGTACCAGGCGGACAGCCCC
CCTGCAGACCTGCTGCTGACGGAGGAGGAGTTCCTGTCGTTCTCCACCCGAGCACAGCCGGGAATGC
TCAGTTTCATGGTGAAGGAGATCGTCCGGACCTGGACCAGGACGGTGACAAGCAGCTCCTGTGCCCGA
GTTTCATCTCCCTGCCCGTGGGCACCGTGGAGAACCAGCAGGGCCAGGACATTGACGACAACCTGGGTGAAA
GACAGAAAAAAGGAGTTTGGAGGAGTCAATTGACTCCAACCACGACGGCATCGTGACCGCCGAGGAGCTGG
AGAACGTGCCCACTCCCGCTGCAGCCAATAGGCACCTAAATAGCCACTTCGTGCGGCTGGCCGCGGA
GCTCGGAGGGGAAAGGCGACGCTGACCTGTGCCCGCTCGCCCGAGAGCTACATGGACCCCA

AG**GCGACCG**ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCC
TGGATTACAAGGATGACGACGATAAGGTTTAA



[View online »](#)


Protein Sequence: >RC216896 representing NM_016547
Red=Cloning site Green=Tags(s)

MVWPWVAMASRWGPLIGLAPCCLLWLLGAVLLMDASARPANHSSTRERVANREENEILPPDHLNGVKLEMD
 GHLNRGFHQEVFLGKDLGGFDEDAEPRRSRRKLMVIFSKVDVNTDRKISAKEMQRWIMEKTAEHFQEAME
 ESKTHFRAVDPDGDGHVSWDEYKVKFLASKGHSEKEVADAIRLNEELKYDEETQEVLENLKDRWYQADSP
 PADLLL TEEEFLSFLHPEHSRGMRLRFMVKEIVRDLDDQDGDQKLSVPEFISLPGTVENQQGQDIDDNWK
 DRKKEFEELIDSNHDGIVTAEELENVPTLPLQPIGTLNSHFVRLAAELGGGKATLTCAPLARRATWTP


SGPTRTRRLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-RsrII

Cloning Scheme:



Plasmid Map:



ACCN: NM_016547

ORF Size: 1044 bp

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_016547.2 , NP_057631.1
RefSeq Size:	2167 bp
RefSeq ORF:	1026 bp
Locus ID:	51150
UniProt ID:	Q9BRK5
Cytogenetics:	1p36.33
Domains:	EFh
Protein Families:	Transmembrane
MW:	39.6 kDa
Gene Summary:	This gene encodes a stromal cell derived factor that is a member of the CREC protein family. The encoded protein contains six EF-hand motifs and calcium-binding motifs. This protein localizes to the Golgi lumen and may be involved in regulating calcium dependent cellular activities. [provided by RefSeq, Sep 2011]