

Product datasheet for **RC216897L3V**

KCNH2 (NM_172057) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	KCNH2 (NM_172057) Human Tagged ORF Clone Lentiviral Particle
Symbol:	KCNH2
Synonyms:	ERG-1; ERG1; H-ERG; HERG; HERG1; Kv11.1; LQT2; SQT1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_172057
ORF Size:	2457 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC216897).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_172057.2 , NP_742054.1
RefSeq Size:	3212 bp
RefSeq ORF:	2460 bp
Locus ID:	3757
UniProt ID:	Q12809
Cytogenetics:	7q36.1
Domains:	cNMP, ion_trans
Protein Families:	Druggable Genome, Ion Channels: Potassium, Transcription Factors, Transmembrane



[View online »](#)

MW: 90.1 kDa

Gene Summary: This gene encodes a voltage-activated potassium channel belonging to the eag family. It shares sequence similarity with the *Drosophila* ether-a-go-go (*eag*) gene. Mutations in this gene can cause long QT syndrome type 2 (LQT2). Transcript variants encoding distinct isoforms have been identified. [provided by RefSeq, Jul 2008]