

Product datasheet for **RC216908A1V**

Human Ovary specific acidic protein (MGARP) (NM_032623) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human Ovary specific acidic protein (MGARP) (NM_032623) AAV Particle
Tag:	Myc-DDK
Symbol:	Ovary specific acidic protein
Synonyms:	C4orf49; CESP-1; HUMMR; OSAP
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC216908 representing NM_032623 Red=Cloning site Blue=ORF Green=Tags(s) TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC ATGTATCTCCGCAGGGCGGTCTCCAAGACTCTGGCGCTGCCGCTGAGGGCGCCCCCAACCCGCGCCGC TCGAAAGGACGCATCTCTGCGCCGGATGTCATCTAACAGATTCCTGGATCATCTGGATCAAATATGAT TTATTATCTGGTTGTAGGCGTCACAGTCAGTGCTGGTGGATATTATGCTTACAAGACAGTCACATCAGAC CAAGCCAAACACACAGAACATAAAACAAATTTGAAAGAAAAACAAAAGCAGAGATACATCCATTTCAAG GTGAAAAGGAGAATGTTGCGGAACTGAGAAAGCAAGTTCAGAAGCCCCAGAAGAACTTATAGTGAAGC TGAGGTGGTAGATGCTGAAGAAAGTCCCAGTGCTACAGTTGTGGTCATAAAAGAGGCATCTGCCTGTCCA GGTCACGTGGAGGCTGCTCCGGAGACCACAGCAGTCAGTGCTGAAACCGGGCCAGAGGTCACAGATGCAG CGGCGAGGGAAACCGGAAGTAAACCTGAAACAACCCAGAGGTTACAAATGCTGCCCTGGATGAAGC TGTCACCATCGATAATGATAAAGATAACAAGAAGCAACCTCTGATGAATATGCTGAAGTAGAAGAA GAAAATTCTCCAGCTGAGTCAGAGTCCTCTGCTGGAGATGATTTACAGGAGGAAGCCAGTGTGGCTCTG AGGCTGCTTCGGCTCAAAGGGCGAATCTCCAGCCGGTAGACATTTAGCAACAAATGCCATAGGGTGTCT GATAAGTGCTTTTGTGTTTTAGTACACTTAGTT ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC216908 representing NM_032623 Red=Cloning site Green=Tags(s) MYLRRRAVSKTLALPLRAPNPAPLGKDALRRMSSNRFPGSSGSNMIYYLVVGVTVSAGGGYAYKTVTSD QAKHTEHKTNLKEKTKAEIHPFQGEKENVAETEKASSEAPEELIVEAEVVDAEESPSATVVVIKEASACP GHVEAAPETTAVSAETGPEVTDAAARETTEVNPETTPETVNAALDEAVTIDNDKDTTKNETSDEYAELEE ENSPAESSAGDDLQEEASVGSEAASAQRANLQPVDISATNAIGCLISAFVFLVHLV TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_032623
ORF Size:	720 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_032623.3</u>
RefSeq Size:	1219 bp
RefSeq ORF:	723 bp
Locus ID:	84709
UniProt ID:	<u>Q8TDB4</u>
Cytogenetics:	4q31.1
MW:	25.2 kDa