

## Product datasheet for RC216917L4V

### STEAP4 (NM\_024636) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	STEAP4 (NM_024636) Human Tagged ORF Clone Lentiviral Particle
Symbol:	STEAP4
Synonyms:	SchLAH; STAMP2; TIARP; TNFAIP9
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_024636
ORF Size:	1377 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC216917).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_024636.2</a>
RefSeq Size:	4454 bp
RefSeq ORF:	1380 bp
Locus ID:	79689
UniProt ID:	<a href="#">Q687X5</a>
Cytogenetics:	7q21.12
Domains:	F420_oxidored
Protein Families:	Druggable Genome, Transmembrane



[View online »](#)

**MW:** 51.8 kDa

**Gene Summary:** The protein encoded by this gene belongs to the STEAP (six transmembrane epithelial antigen of prostate) family, and resides in the golgi apparatus. It functions as a metalloreductase that has the ability to reduce both Fe(3+) to Fe(2+) and Cu(2+) to Cu(1+), using NAD(+) as acceptor. Studies in mice and human suggest that this gene maybe involved in adipocyte development and metabolism, and may contribute to the normal biology of the prostate cell, as well as prostate cancer progression. Alternatively spliced transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Apr 2011]