

Product datasheet for **RC216931L3V**

CDC14A (NM_033312) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	CDC14A (NM_033312) Human Tagged ORF Clone Lentiviral Particle
Symbol:	CDC14A
Synonyms:	cdc14; DFNB32; DFNB35; DFNB105; hCDC14
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_033312
ORF Size:	1869 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC216931).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_033312.1 , NP_201569.1
RefSeq Size:	2438 bp
RefSeq ORF:	1872 bp
Locus ID:	8556
UniProt ID:	Q9UNH5
Cytogenetics:	1p21.2
Domains:	Y_phosphatase
Protein Families:	Druggable Genome, Phosphatase



[View online »](#)

Protein Pathways: Cell cycle

MW: 69.4 kDa

Gene Summary: The protein encoded by this gene is a member of the dual specificity protein tyrosine phosphatase family. It is highly similar to *Saccharomyces cerevisiae* Cdc14, a protein tyrosine phosphatase involved in the exit of cell mitosis and initiation of DNA replication, suggesting a role in cell cycle control. This protein has been shown to interact with, and dephosphorylate tumor suppressor protein p53, and is thought to regulate the function of p53. Alternative splicing of this gene results in several transcript variants encoding distinct isoforms. [provided by RefSeq, Jul 2008]