

## Product datasheet for RC216947A1V

### Human IL13 (NM\_002188) AAV Particle

#### Product data:

**Product Type:** AAV Particles  
**Product Name:** Human IL13 (NM\_002188) AAV Particle  
**Tag:** Myc-DDK  
**Symbol:** IL13  
**Synonyms:** IL-13; P600  
**Mammalian Cell Selection:** None  
**Vector:** pAAV-AC-Myc-DDK (PS100089)  
**ORF Nucleotide Sequence:** >RC216947 representing NM\_002188  
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGGATCGCC**

ATGCATCCGCTCCTCAATCCTCTCCTGTTGGCACTGGCCTCATGGCGCTTTTGTGACCACGGTCATTG  
CTCTCACTTGCCTTGGCGGCTTTGCCTCCCCAGGCCCTGTGCCTCCCTCTACAGCCCTCAGGGAGCTCAT  
TGAGGAGCTGGTCAACATCACCCAGAACCAGAAGGCTCCGCTCTGCAATGGCAGCATGGTATGGAGCATC  
AACCTGACAGCTGGCATGTACTGTGCAGCCCTGGAATCCCTGATCAACGTGTCAGGCTGCAGTGCCATCG  
AGAAGACCCAGAGGATGCTGAGCGGATTCTGCCCGCACAAAGGTCTCAGCTGGGCAGTTTTCCAGCTTGCA  
TGTCCGAGACACCAAAATCGAGGTGGCCAGTTTGTAAGGACCTGCTTACATTTAAAGAAACTTTTT  
CGCGAGGGACAGTTCAAC

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:** >RC216947 representing NM\_002188  
 Red=Cloning site Green=Tags(s)

MHPLLNP LLLALGLMALLLTTVIALTCLGGFASPGVPVPPSTALRELIEELVNITQNQKAPLCNGSMVWSI  
NLTAGMYCAALESLINVSGCSAIEKTQRMLSGFCPHKVSAGQFSSLHVRDTKIEVAQFVKDLLLHLKLF  
REGQFN

**TR**TRPLEQKLI SEEDLAANDILDYKDDDDKV

**Species:** Human  
**Serotype:** AAV-2



[View online »](#)

ACCN:	NM_002188
ORF Size:	438 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u><a href="#">NM_002188.2</a></u>
RefSeq Size:	1282 bp
RefSeq ORF:	441 bp
Locus ID:	3596
UniProt ID:	<u><a href="#">P35225</a></u>
Cytogenetics:	5q31.1
MW:	15.79 kDa