

## Product datasheet for RC216953

### OR10J1 (NM\_012351) Human Tagged ORF Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** OR10J1 (NM\_012351) Human Tagged ORF Clone  
**Tag:** Myc-DDK  
**Symbol:** OR10J1  
**Synonyms:** HGMP07J; HSHGMP07J  
**Mammalian Cell Selection:** Neomycin  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**ORF Nucleotide Sequence:** >RC216953 representing NM\_012351  
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGCATCGCC**

ATGCTTTTATGTTTCAGATTTGGGAACCAATCCATGAAAAGAGAGAACTTTACTCTCATCACTGACTTTG  
TTTTCCAAGGTTTCTCTAGCTTCCATGAGCAGCAGATCACCCTTTTGGCGTGTTCTTGCACTATAACAT  
CTTAACCTTAGCAGGCAATATCATCATTGTGACCATCATCCGAATGGATTTTCATCTTACACACCCATG  
TACTTCTTCTGAGCATGCTGTCCACTTCAGAGACTGTATACATTGGTCATTCTCCAAGAATGCTCT  
CCAGCCTCGTAGGTATGAGCCAGCCCATATCATTGGCAGGGTGTCCACACAGATGTTCTTTTTGTAAC  
CTTTGGCATCACTAACTGCTTCCGTCTCACAGCAATGGGATATGACCGCTATGTGGCCATCTGCAACCC  
CTGAGATACATGGTTATTATGAACAAGAGGCTGCGTATCCAACCTGTCTGGGGGCTGCAGCATTGGGC  
TGATTGTAGCAATAACGCAAGTGACATCTGTATTACAGTTACCCTTCTGTGCTAGAAAGGTGCCCACTT  
CTTCTGTGACATCCGCCCTGTGATGAAGCTCTCCTGCAATTGACACCACTGTCAATGAAATCCTGACTTTG  
ATTATCAGTGTGCTGGTCTTGTGTACCTATGGGTCTGGTTTTTCATTTCTTATGTTCTCATTATCTCTA  
CAATCCTCAAGATTGCTTCAGTTGAGGGCCGGAAGAAGCTTTTGCCACCTGTGCATCCCACCTCACTGT  
GGTCATTGTCCACTACAGCTGTGCCTCCATTGCCTACCTCAAGCCCAAGTCAGAGAACACCAGAGAACAT  
GACCAGCTGATCTCGGTGACCTACACTGTATCACTCCCCTACTGAACCCCTGTGGTATACACCCTGAGAA  
ATAAAGAGGTCAAAGATGCTCTGTGCAGGGCTGTTGGTGGGAAGTTTTCC

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

**Protein Sequence:** >RC216953 representing NM\_012351  
Red=Cloning site Green=Tags(s)

MLLCFRFGNQSMKRENFTLITDFVFQGFSSFHEQQITLFGVFLALYLTLAGNIIIVTIIRMDLHLHTPM  
 YFFLSMLSTSETVYTLVILPRMLSSLVGMSQPISLAGCATQMFVFTFGITNCFLLTAMGYDRYVAICNP  
 LRYMVMINKRLRIQLVLGACISGLIVAITQVTSVFRLPFCARKVPHFCDIRPVMKLSCIDTTVNEILTL  
 IISVLVLVPMGLVFIISVLIISTILKIASVEGRKKAFATCASHLTVVIVHYSCASIAYLKPKSENTREH  
 DQLISVTTYVITPLLNPVYVYTLRNKEVKDALCRAVGGKFS

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

**Chromatograms:** [https://cdn.origene.com/chromatograms/mk8008\\_g02.zip](https://cdn.origene.com/chromatograms/mk8008_g02.zip)

**Restriction Sites:** SgfI-MluI

**Cloning Scheme:**

Cloning sites used for ORF Shuttling:



\* The last codon before the Stop codon of the ORF

**ACCN:** NM\_012351

**ORF Size:** 960 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_012351.2](#), [NP\\_036483.2](#)

**RefSeq Size:** 963 bp

**RefSeq ORF:** 930 bp

**Locus ID:** 26476

**UniProt ID:** [P30954](#)

**Cytogenetics:** 1q23.2


**Protein Families:** Druggable Genome, Transmembrane

**Protein Pathways:** Olfactory transduction

**MW:** 35.7 kDa

**Gene Summary:** Olfactory receptors interact with odorant molecules in the nose, to initiate a neuronal response that triggers the perception of a smell. The olfactory receptor proteins are members of a large family of G-protein-coupled receptors (GPCR) arising from single coding-exon genes. Olfactory receptors share a 7-transmembrane domain structure with many neurotransmitter and hormone receptors and are responsible for the recognition and G protein-mediated transduction of odorant signals. The olfactory receptor gene family is the largest in the genome. The nomenclature assigned to the olfactory receptor genes and proteins for this organism is independent of other organisms. [provided by RefSeq, Jul 2008]

Product images:



Circular map for RC216953