

Product datasheet for **RC216956L4V**

Deoxyribonuclease I like 1 (DNASE1L1) (NM_001009932) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Deoxyribonuclease I like 1 (DNASE1L1) (NM_001009932) Human Tagged ORF Clone Lentiviral Particle
Symbol:	Deoxyribonuclease I like 1
Synonyms:	DNAS1L1; DNASEX; DNL1L; G4.8; XIB
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001009932
ORF Size:	906 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC216956).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001009932.1 , NP_001009932.1
RefSeq Size:	3008 bp
RefSeq ORF:	909 bp
Locus ID:	1774
UniProt ID:	P49184
Cytogenetics:	Xq28
MW:	33.9 kDa



[View online »](#)

Gene Summary:

This gene encodes a deoxyribonuclease protein that shows high sequence similarity to DNase I. The encoded protein is localized to the endoplasmic reticulum and modified by N-linked glycosylation. Alternate transcriptional splice variants encoding the same protein have been observed. [provided by RefSeq, Jan 2015]