

Product datasheet for **RC216958A1V**

Human JUND (NM_005354) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human JUND (NM_005354) AAV Particle
Tag:	Myc-DDK
Symbol:	JUND
Synonyms:	AP-1
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC216958 representing NM_005354 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGGAACACCCCTTCTACGGCGATGAGGCGCTGAGCGGCTGGGCGGCGGCCAGTGCCAGCGGCGGCA
GCTTCGCGTCCCCGGGCGCTTGTTCGCCGGGCGCCCCGACGGCCGCGGCCGCGCAGCATGATGAAGAA
GGACGCGCTGACGCTGAGCCTGAGTGAGCAGGTGGCGGCAGCGCTCAAGCCTGCGGCCGCGCCCTCCT
ACCCCCCTGCGCGCCGACGGCGCCCCAGCGCGGCACCCCCGACGGCCTGCTCGCCTCTCCCGACCTGG
GGCTGCTGAAGCTGGCTCCCCGAGCTCGAGCGCCTCATCATCCAGTCCAACGGGCTGGTCACCAACCAC
GCCGACGAGCTCACAGTTCTCTACCCCAAGGTGGCGGCCAGCGAGGAGCAGGAGTTCCGCCGAGGGCTTC
GTCAAGGCCCTGGAGGATTTACACAAGCAGAACCAGCTCGGCGCGGGCGCGGCCGCTGCCGCCGCGCCG
CCGCCGCCGGGGGCCCTCGGGCACGGCCACGGGCTCCGCGCCCCCGGCGAGCTGGCCCCGGCGCGGC
CGCGCCGAAGCGCCTGTCTACGCGAACCTGAGCAGCTACGCGGGCGGCGCCGGGGGCGCGGGGGCGCC
GCGACGGTCGCTTCTGCTGCCGAACCTGTGCCCTTCCGCGCCGCCACCCCCAGGCGCGTTGGGGCCGC
CGCGCCTGGCTGCGCTCAAGGACGAGCCACAGACGGTGCCCGACGTGCCGAGCTTCGGCGAGAGCCCGCC
GTTGTCGCCCATCGACATGGACACGCAGGAGCGCATCAAGGCGGAGCGCAAGCGGCTGCGCAACCGCATC
GCCGCCTCCAAGTGCCGAAGCGCAAGCTGGAGCGCATCTCGCGCCTGGAAGAGAAAGTGAAGACCTCA
AGAGTCAGAACCGGAGCTGGCGTCCACGGCGAGCCTGCTGCGCGAGCAGGTGGCGCAGCTCAAGCAGAA
AGTCTCAGCCAGCTAACAGCGGCTGCCAGCTGCTGCCCCAGCACCAGGTGCCCGCGTAC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC216958 representing NM_005354 Red=Cloning site Green=Tags(s) METPFYGDEALSGLGGGASGSGGSFASPGRLFPGAPPTAAAGSMMKKDALTLSEQVAAALKPAAAPPP TPLRADGAPSAAPPDGLLASPDGLLLKASPELERLIIQSNGLVTTTPTSSQFLYPKVAASEEQEFAEGF VKALEDLHKQNLGAGAAAAAAAAAGGPGSGTATGSAPPGELAPAAAAPEAPVYANLSSYAGGAGGAGGA ATVAFAAEPVPFPPPPPGALGPPRLAALKDEPQTVPDVPSFGESPLSPIDMDTQERIKAEKRLRNRI AASKCRKRKLERISRLEEKVKTLKSQNTELASTALLREQVAQLKQKVLSHVNSGCQLLPQHQPAY TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_005354
ORF Size:	1041 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_005354.2
RefSeq Size:	1891 bp
RefSeq ORF:	1044 bp
Locus ID:	3727
UniProt ID:	P17535
Cytogenetics:	19p13.11
MW:	35 kDa