

Product datasheet for **RC216982A1V**

Human CD62L (SELL) (NM_000655) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human CD62L (SELL) (NM_000655) AAV Particle
Tag:	Myc-DDK
Symbol:	CD62L
Synonyms:	CD62L; LAM1; LECAM1; LEU8; LNHR; LSEL; LYAM1; PLNHR; TQ1
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC216982 representing NM_000655 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGATATTTCCATGGAATGTCAGAGCACCCAGAGGGACTTATGGAACATCTTCAAGTTGTGGGGTGGAC
CAATGCTCTGTTGTGATTTCTGGCACATCATGGAACCGACTGCTGGACTTACCATTATTTCTGAAAAACC
CATGAAGTGGCAAAGGGCTAGAAGATTCTGCCGAGACAATTACACAGATTTAGTTGCCATACAAAAACAAG
GCGGAAATTGAGTATCTGGAGAAGACTCTGCCTTCAGTCGTTCTTACTACTGGATAGGAATCCGGAAGA
TAGGAGGAATATGGAGCTGGGTGGGAACCAACAATCTTACTGAAGAAGCAGAGAAGTGGGGAGATGG
TGAGCCCAACAACAAGAAGAACAAGGAGGACTGCGTGGAGATCTATATCAAGAGAAACAAGATGCAGGC
AAATGGAACGATGACGCCTGCCACAACTAAAGGCAGCCCTCTGTTACACAGCTTCTTGCCAGCCCTGGT
CATGCAGTGGCCATGGAGAATGTGTAGAAATCATCAATAATTACACCTGCAACTGTGATGTGGGGTACTA
TGGGCCCCAGTGTGAGTTTGTGATTCAGTGTGAGCCTTTGGAGGCCCCAGAGCTGGGTACCATGGACTGT
ACTCACCTTTGGGAACTTCAGCTTCAGCTCACAGTGTGCCTTCAGCTGCTCTGAAGGAACAACTTAA
CTGGGATTGAAGAAACCACTGTGGACCATTTGGAACTGGTCATCTCCAGAACCAACCTGTCAAGTGT
TCAGTGTGAGCCTCTATCAGCACCAGATTTGGGGATCATGAAGTGTAGCCATCCCCTGGCCAGCTTCAGC
TTTACCTCTGCATGTACCTTCATCTGCTCAGAAGGAAGTGAATTAATTGGGAAGAAGAAACCATTTGTG
AATCATCTGGAATCTGGTCAAACTCTAGTCCAATATGTCAAAAATTGGACAAAAGTTTCTCAATGATTAA
GGAGGGTGATTATAACCCCTCTTCATTCCAGTGGCAGTCATGTTTACTGCATTCTCTGGTTGGCATT
ATCATTGGCTGGCAAGGAGATTAATAAAGGCAAGAAATCCAAGAGAAGTATGAATGACCCATAT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC216982 representing NM_000655 Red=Cloning site Green=Tags(s) MIFPWKCQSTQRDLWNIFKLWGWTMLCCDFLAHHGTDWCWYHYSEKPMNWQRARRFCRDNYTDLVAIQNK AEIEYLEKTLPFSSRSYYWIGIRKIGGIWTWVGTNKSLEEAENWGDGEPNNKKNKEDCVEIYIKRNKDAG KWND DACHKLKAALCYTASCQPWSCSGHGECVEIINNYTCNCDVGYGPGQCQFVIQCEPLEAPELGTMD THPLGNFSFSSQCAFSCSEGTNLTGIEETTCGPFGNWSSPEPTCQVIQCEPLSAPDLGIMNC SHPLASFS FTSACTFICSEGTELIGKKKTICSSGIWSNPSPICQKLDKSF SMIKEGDYNPLFIPVAVMVTAFSGLAF IIWLARRLKKGKSKRSMNDPY TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_000655
ORF Size:	1116 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_000655.2
RefSeq Size:	2324 bp
RefSeq ORF:	1119 bp
Locus ID:	6402
UniProt ID:	P14151
Cytogenetics:	1q24.2
MW:	42.19 kDa