

## Product datasheet for **RC216997L3V**

### **BTN2A1 (NM\_007049) Human Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	BTN2A1 (NM_007049) Human Tagged ORF Clone Lentiviral Particle
Symbol:	BTN2A1
Synonyms:	BK14H9.1; BT2.1; BTF1; BTN2.1; DJ3E1.1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_007049
ORF Size:	1581 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC216997).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_007049.3</a>
RefSeq Size:	2906 bp
RefSeq ORF:	1584 bp
Locus ID:	11120
UniProt ID:	<a href="#">Q7KYR7</a>
Cytogenetics:	6p22.2
Domains:	IGv, IG, SPRY, PRY
Protein Families:	Druggable Genome, Transmembrane




[View online »](#)


MW: 59.63 kDa

**Gene Summary:** This gene encodes a member of the immunoglobulin superfamily. The gene is located in a cluster of butyrophilin-like genes in the juxta-telomeric region of the major histocompatibility complex on chromosome 6. A pseudogene of this gene has been identified in this cluster. The encoded protein is an integral plasma membrane protein involved in lipid, fatty-acid, and sterol metabolism. Alterations in this gene may be associated with several disease states including metabolic syndrome. Multiple alternatively spliced transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jul 2013]

### Product images:



[RC216997L3] was used to prepare Lentiviral particles using [TR30037] packaging kit. HEK293T cells were transduced with RC216997L3V particle to overexpress human BTN2A1-Myc-DDK fusion protein.



[RC216997L3] was used to prepare Lentiviral particles using [TR30037] packaging kit. HEK293T cells were transduced with RC216997L3V particle to overexpress human BTN2A1-Myc-DDK fusion protein.