

## Product datasheet for RC217032L4V

### CRCP (NM\_014478) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	CRCP (NM_014478) Human Tagged ORF Clone Lentiviral Particle
Symbol:	CRCP
Synonyms:	C17; CGRP-RCP; CGRPRCP; POLR3I; POLR3J; RCP; RCP9; RPC9
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_014478
ORF Size:	444 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC217032).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_014478.3</a>
RefSeq Size:	1048 bp
RefSeq ORF:	447 bp
Locus ID:	27297
UniProt ID:	<a href="#">O75575</a>
Cytogenetics:	7q11.21
Domains:	RPOL4c
Protein Families:	Druggable Genome



[View online »](#)

**MW:** 16.7 kDa

**Gene Summary:** This gene encodes a membrane protein that functions as part of a receptor complex for a small neuropeptide that increases intracellular cAMP levels. Alternate transcriptional splice variants, encoding different isoforms, have been characterized. [provided by RefSeq, Jul 2008]