

Product datasheet for **RC217042A1V**

Human LINC02913 (NM_001012715) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Human LINC02913 (NM_001012715) AAV Particle
Tag: Myc-DDK
Symbol: LINC02913
Synonyms: bA65J3.5
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC217042 ORF sequence
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGGCGTATGTGGCCTTGAGCGATAAACCTCACCTGTCGGGGGAGGTGGGCGAGGAAGCCTGCAGCTCCT
GGAACCCACCACTTTTCTCTCCAGGCCCTTGCCAGGCTGTGGCCTCTGCCTGGAACCCCTTTTCTTCC
CATCCGAGCCCTTCCCTTCAGTGCCCTCTTCTCTGGGAAGTCTCTCTGGTGGCCAGCCCGTCTCTTCT
GCCACTCTGGGTTTCGCACCCCATGCTTAGGCCCTGATTGCCTGCTTTGCACACAGGGCTGTGAGCTCC
ATGAGGGCAGGAACCATGGCTGTTTACAGCTGTGTTGCCAGGGCCTGGCCTGGAGACCCTCAGGAAGT
GAGACATCTGAATCCTCTGCTCTGTGACCCCTGGCAGCCAGGTAGAGCCATCTTGCCCATGGCAGCCAGGC
TTGGAGCAAGCAGCAGCCTCATGGGTGGGAACACGTTTCCCAGCACATCGGCAGGCCTAAGGGGAC
ATTCGCTTGGCTCTGCCCTCAGAGCACTGATGCCCGGAGACTGCCCTTTGTGTGCCATGTAAGAG
AGGATGCAACCTTAGAGGGGCGAGGGTAAACATGGTCCAGGCCCTGTTGTCCCTCTAAGAAAGTTC
CCAGTTCTGCCTGTCCATCCCTGGCCTTTTCCCTGTGCTGTCTGGGACAGTGGGTGGAGCTGCCGC

ACGCGTACGCGGCGGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC217042 protein sequence
Red=Cloning site Green=Tags(s)

MAYVALSDKPHLSGEVGEEACSSWNPPFLSPRPLPRLWPLPGTFPLPIRALPFSASSSGKSSLVSPSSS
AHSGRFTPLGPDCLLTQGCHEGRNHMAVHSCVARAWPGDPQEVRLNPLLCDPGSQVEPSWPWHPG
LEQAAASWVGNHVSPAHRQALRGHSLGSALRALMPGRHCPLCVPCKRGCNLRGGRGKHGPRPCCPLLRF
PVLPHWPFPFCAVWDSGWSCR

TRTRPLEQKLISEEDLAANDILDYKDDDDKV



[View online »](#)

Species:	Human
Serotype:	AAV-2
ACCN:	NM_001012715
ORF Size:	696 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_001012715.1</u> , <u>NP_001012733.2</u>
RefSeq Size:	843 bp
RefSeq ORF:	699 bp
Locus ID:	414318
Cytogenetics:	9q34.11
MW:	25.1 kDa