

Product datasheet for **RC217063A1V**

Human HAUS8 (NM_033417) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human HAUS8 (NM_033417) AAV Particle
Tag:	Myc-DDK
Symbol:	HAUS8
Synonyms:	DGT4; HICE1; NY-SAR-48
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC217063 representing NM_033417 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGGCGGATTCTCGGGCGAGGCGCTGGGAAGCCTGCAACCGGCCCCACAAATTCTAGCAGTGCCAAGA
AGAAGGATAAAAGAGTTCAAGGTGGAAGAGTGATTGAGTCCCGGTATCTGCAGTATGAAAAGAAGACAAC
CCAAAAGGCTCCTGCAGGAGATGGGTACAGACCCGAGGAAGATGTCTGAAGGTGGAAGGAAATCCAGC
CTGCTCCAGAAAAGCAAAGCAGATAGCAGTGGGTGCGAAAGGGTGACCTGCAGTCCACGTTGCTGGAAG
GGCATGGCACAGCTCCACCTGACCTGGATCTCTGCTATTAATGACAAAAGCATCGTCAAAAAGACGCC
ACAGTTAGCAAAAACAATATCAAAAGAAACCTGAGTCAACATCATTTTCTGCCCTCGGAAAAAGAGCCCG
GATTTATCTGAAGCAATGGAATGATGGAGTCTCAGACACTACTGCTGACGCTACTATCCGTAAAGATGG
AGAACAATCTTGCTGAGTTTGAAAGAAGGGCAGAAAAGAATTTATTAATAATGTGTAAGGAGAAGGAGAA
GCTACAGAAAAGGCCACGAGCTGAAGCGCAGGCTTCTCTCTCTCAGAGGAAGCGGGAGCTGGCAGAT
GTCCTGGATGCCAGATCGAGATGCTCAGCCCCCTCGAGGCAGTGGCCACACGCTTCAAGGAGCAATACA
GGACATTCGCCACGGCCCTGGACACTACCAGGCACGAGCTGCCCGTGAGGTCCATCCACCTGGAGGGAGA
TGGGCAGCAGCTCTTAGACGCCCTGCAGCATGAACTGGTGACCACTCAGCGCCTCTGGGAGAACTTGAT
GTTGGTGATTGGAAGAAAATGTGCAGGTGCTGGACTTACTGAGCGAACTCAAGGACGTGACGGCGAAAA
AGGACCTTGAGCTCCGAAGGAGCTTTGCCAGGTGCTGGAACTCTCCGCAGAGGCAAGCAAGAGGCAGC
CTTGGCAAACAGGAAGTCTGGGAAGAGACCCAGGGCATGGCGCCCCCAGCCGGTGGTATTTCATCAA
GACAGTGCCTGCAGAGAATCTGGGGGAGCACCAAGAACACGCCCTGTCTGAGGACGACAACCCGGGTG
CCTCGTCAGCCCCGCTCAGGCCACGTTTCATCAGCCCAAGCGAAGATTTTCTTCAAGCAGCCAGGCAGA
AGTCCCGCCCTCTCTCTCTCTCGTTCAGGGAGGGACTTGTC

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC217063 representing NM_033417 Red=Cloning site Green=Tags(s) MADSSGRGAGKPATGPTNSSSAKKKDKRVQGGRVIESRYLQYEKKTQKAPAGDGSQTRGKMSEGGKSS LLQKSKADSSGVGKDLQSTLLEGHGTAPPDLDSAINDKSIVKKTQQLAKTISKKPESTSFSAPRKKSP DLSEAMEMMESQTLLLTLLSVKMENNLAEFERRAEKNLLIMCKEKEKLQKKAHELKRRLLSQRKRELAD VLDAQIEMLSPF EAVATRFKEQYRTFATALDTTRHEL PVR SIHLEGGQQLLDALQHEL VTTQRLLGELD VGDSEENVQVLDLLSELKDVTAKKDLELRRSFAQVLELSAEASKEAALANQEVWEETQGMAPP SRWYFNQ DSACRESGGAPKNTPLSEDDNPGASSAPAQATFISPS EDFSSSSQAEVPPSLSRSGRDLS TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_033417
ORF Size:	1230 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_033417.1
RefSeq Size:	1468 bp
RefSeq ORF:	1233 bp
Locus ID:	93323
UniProt ID:	Q9BT25
Cytogenetics:	19p13.11
MW:	44.7 kDa