

## Product datasheet for **RC217083L1V**

### SKAP55 (SKAP1) (NM\_003726) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	SKAP55 (SKAP1) (NM_003726) Human Tagged ORF Clone Lentiviral Particle
Symbol:	SKAP55
Synonyms:	HEL-S-81p; SCAP1; SKAP55
Mammalian Cell Selection:	None
Vector:	pLenti-C-Myc-DDK (PS100064)
Tag:	Myc-DDK
ACCN:	NM_003726
ORF Size:	1077 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC217083).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_003726.2</a> , <a href="#">NP_003717.2</a>
RefSeq Size:	1531 bp
RefSeq ORF:	1080 bp
Locus ID:	8631
UniProt ID:	<a href="#">Q86WV1</a>
Cytogenetics:	17q21.32
Domains:	SH3, PH
MW:	41.3 kDa


[View online »](#)

**Gene Summary:**

This gene encodes a T cell adaptor protein, a class of intracellular molecules with modular domains capable of recruiting additional proteins but that exhibit no intrinsic enzymatic activity. The encoded protein contains a unique N-terminal region followed by a PH domain and C-terminal SH3 domain. Along with the adhesion and degranulation-promoting adaptor protein, the encoded protein plays a critical role in inside-out signaling by coupling T-cell antigen receptor stimulation to the activation of integrins. [provided by RefSeq, Jul 2008]