

Product datasheet for **RC217096L3V**

zinc finger protein 655 (ZNF655) (NM_001085366) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	zinc finger protein 655 (ZNF655) (NM_001085366) Human Tagged ORF Clone Lentiviral Particle
Symbol:	ZNF655
Synonyms:	VIK; VIK-1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001085366
ORF Size:	330 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC217096).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001085366.1 , NP_001078835.1
RefSeq Size:	1395 bp
RefSeq ORF:	333 bp
Locus ID:	79027
UniProt ID:	Q8N720
Cytogenetics:	7q22.1
Protein Families:	Transcription Factors



[View online »](#)

MW: 11.9 kDa

Gene Summary: This gene encodes a zinc finger protein. The zinc finger proteins are involved in DNA binding and protein-protein interactions. Multiple alternatively spliced transcript variants encoding distinct isoforms have been found for this gene. [provided by RefSeq, Jul 2008]