

## Product datasheet for **RC217103**

### Glycophorin C (GYPC) (NM\_016815) Human Tagged ORF Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Glycophorin C (GYPC) (NM\_016815) Human Tagged ORF Clone  
**Tag:** Myc-DDK  
**Symbol:** Glycophorin C  
**Synonyms:** CD236; CD236R; GE; GE:GPC:GPD:GYPD; GPC; GPD; GYPD; PAS-2; PAS-2'  
**Mammalian Cell Selection:** Neomycin  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**ORF Nucleotide Sequence:** >RC217103 representing NM\_016815  
**Red=Cloning site Blue=ORF Green=Tags(s)**

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGTGGTCGACGAGAAGCCCCAACAGCACGGCGTGGCCTCTCAGCCTCGAGCCTGATCCAGGGATGTCTG  
GATGGCCGGATGGCAGAATGGAGACCTCCACCCACCATAATGGACATTGTCGTCATTGCAGGTGTGAT  
TGCTGCTGTGGCCATCGTCCTAGTCTCCCTCCTCTTCGTCATGCTGCGCTACATGTACCGGCACAAGGGC  
ACGTACCACCAATGAGGCCAAGGGCACGGAGTTTGCTGAGAGTGCAGATGCAGCCCTGCAGGGAGACC  
CTGCCCTCAAGATGCTGGTATAGCAGCAGAAAGGAGTACTTTATT

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:** >RC217103 representing NM\_016815  
**Red=Cloning site Green=Tags(s)**  
MWSTRSPNSTAWPLSLEPDPMGSGWPDGRMETSTPTIMDIVVIAGVIAAVAIVLVSLLFVMLRYMYRHKG  
TYHTNEAKGTEFAESADAALQGDPALQDAGDSSRKEYFI

**TR**TRPLEQKLISEEDLAANDILDYKDDDDKV

**Chromatograms:** [https://cdn.origene.com/chromatograms/mg3409\\_g01.zip](https://cdn.origene.com/chromatograms/mg3409_g01.zip)

**Restriction Sites:** SgfI-MluI



[View online »](#)

**Cloning Scheme:**


**ACCN:** NM\_016815

**ORF Size:** 327 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_016815.3](#)

**RefSeq Size:** 1019 bp


**RefSeq ORF:** 330 bp

**Locus ID:** 2995


**UniProt ID:** [P04921](#)  
**Cytogenetics:** 2q14.3  
**Protein Families:** Druggable Genome, Transmembrane  
**MW:** 11.7 kDa

**Gene Summary:** Glycophorin C (GYPC) is an integral membrane glycoprotein. It is a minor species carried by human erythrocytes, but plays an important role in regulating the mechanical stability of red cells. A number of glycophorin C mutations have been described. The Gerbich and Yus phenotypes are due to deletion of exon 3 and 2, respectively. The Webb and Duch antigens, also known as glycophorin D, result from single point mutations of the glycophorin C gene. The glycophorin C protein has very little homology with glycophorins A and B. Alternate splicing results in multiple transcript variants. [provided by RefSeq, Feb 2012]


### Product images:



Circular map for RC217103



Western blot validation of overexpression lysate (Cat# [LY413831]) using anti-DDK antibody (Cat# [TA50011-100]). Left: Cell lysates from untransfected HEK293T cells; Right: Cell lysates from HEK293T cells transfected with RC217103 using transfection reagent MegaTran 2.0 (Cat# [TT210002]).



Coomassie blue staining of purified GYPC protein (Cat# [TP317103]). The protein was produced from HEK293T cells transfected with GYPC cDNA clone (Cat# RC217103) using MegaTran 2.0 (Cat# [TT210002]).