

Product datasheet for **RC217115L3V**

MEK5 (MAP2K5) (NM_002757) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	MEK5 (MAP2K5) (NM_002757) Human Tagged ORF Clone Lentiviral Particle
Symbol:	MEK5
Synonyms:	HsT17454; MAPKK5; MEK5; PRKMK5
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_002757
ORF Size:	1314 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC217115).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_002757.2
RefSeq Size:	2334 bp
RefSeq ORF:	1317 bp
Locus ID:	5607
UniProt ID:	Q13163
Cytogenetics:	15q23
Domains:	PB1, pkinase, TyrKc, S_TKc
Protein Families:	Druggable Genome, Protein Kinase



[View online »](#)

Protein Pathways: Gap junction, MAPK signaling pathway, Neurotrophin signaling pathway

MW: 48.8 kDa

Gene Summary: The protein encoded by this gene is a dual specificity protein kinase that belongs to the MAP kinase kinase family. This kinase specifically interacts with and activates MAPK7/ERK5. This kinase itself can be phosphorylated and activated by MAP3K3/MEKK3, as well as by atypical protein kinase C isoforms (aPKCs). The signal cascade mediated by this kinase is involved in growth factor stimulated cell proliferation and muscle cell differentiation. Three alternatively spliced transcript variants of this gene encoding distinct isoforms have been described. [provided by RefSeq, May 2011]