

Product datasheet for RC217139L3

ZNF365 (NM_199450) Human Tagged Lenti ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ZNF365 (NM_199450) Human Tagged Lenti ORF Clone
Tag:	Myc-DDK
Symbol:	ZNF365
Synonyms:	Su48; UAN; ZNF365D
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
E. coli Selection:	Chloramphenicol (34 ug/mL)
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC217139).
Restriction Sites:	SgfI-MluI
Cloning Scheme:	

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF.


ACCN:	NM_199450
ORF Size:	999 bp



[View online »](#)

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_199450.1 , NP_955522.1
RefSeq Size:	1671 bp
RefSeq ORF:	1002 bp
Locus ID:	22891
UniProt ID:	Q70YC5
Cytogenetics:	10q21.2
MW:	38 kDa
Gene Summary:	This gene encodes a zinc finger protein that may play a role in the repair of DNA damage and maintenance of genome stability. The N-terminal C2H2 zinc finger motif is required to form a protein complex with PARP1 and MRE11, which are known to be involved in the restart of stalled DNA replication forks. A mutation in this gene may be associated with breast cancer susceptibility. [provided by RefSeq, Mar 2020]

Product images:



Circular map for RC217139L3