

Product datasheet for **RC217163L3V**

MFGE8 (NM_005928) Human Tagged ORF Clone Lentiviral Particle

Product data:


Product Type:	Lentiviral Particles
Product Name:	MFGE8 (NM_005928) Human Tagged ORF Clone Lentiviral Particle
Symbol:	MFGE8
Synonyms:	BA46; EDIL1; HMFG; hP47; HsT19888; MFG-E8; MFGM; OAcGD3S; SED1; SPAG10
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_005928
ORF Size:	1161 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC217163).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_005928.1
RefSeq Size:	1934 bp
RefSeq ORF:	1164 bp
Locus ID:	4240
UniProt ID:	Q08431
Cytogenetics:	15q26.1
Domains:	F5_F8_type_C, EGF, EGF
MW:	42.9 kDa



[View online »](#)

Gene Summary:

This gene encodes a preproprotein that is proteolytically processed to form multiple protein products. The major encoded protein product, lactadherin, is a membrane glycoprotein that promotes phagocytosis of apoptotic cells. This protein has also been implicated in wound healing, autoimmune disease, and cancer. Lactadherin can be further processed to form a smaller cleavage product, medin, which comprises the major protein component of aortic medial amyloid (AMA). Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jul 2015]

Product images:

[RC217163L3] was used to prepare Lentiviral particles using [TR30037] packaging kit. HEK293T cells were transduced with RC217163L3V particle to overexpress human MFGE8-Myc-DDK fusion protein.