

## Product datasheet for **RC217173L4V**

### **FXR1 (NM\_005087) Human Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	FXR1 (NM_005087) Human Tagged ORF Clone Lentiviral Particle
Symbol:	FXR1
Synonyms:	FXR1P; MYOPMIL; MYORIBF
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_005087
ORF Size:	1863 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC217173).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_005087.2</a>
RefSeq Size:	2774 bp
RefSeq ORF:	1866 bp
Locus ID:	8087
UniProt ID:	<a href="#">P51114</a>
Cytogenetics:	3q26.33
Domains:	KH
Protein Families:	Druggable Genome



[View online »](#)

**MW:** 69.5 kDa

**Gene Summary:** The protein encoded by this gene is an RNA binding protein that interacts with the functionally-similar proteins FMR1 and FXR2. These proteins shuttle between the nucleus and cytoplasm and associate with polyribosomes, predominantly with the 60S ribosomal subunit. Three transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jul 2008]