

## Product datasheet for **RC217183L1V**

### **RNF13 (NM\_007282) Human Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	RNF13 (NM_007282) Human Tagged ORF Clone Lentiviral Particle
Symbol:	RNF13
Synonyms:	DEE73; EIEE73; RZF
Mammalian Cell Selection:	None
Vector:	pLenti-C-Myc-DDK (PS100064)
Tag:	Myc-DDK
ACCN:	NM_007282
ORF Size:	1143 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC217183).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_007282.3</a> , <a href="#">NP_009213.1</a>
RefSeq Size:	2956 bp
RefSeq ORF:	1146 bp
Locus ID:	11342
UniProt ID:	<a href="#">O43567</a>
Cytogenetics:	3q25.1
Domains:	RING, PA
Protein Families:	Druggable Genome, Protease, Transmembrane



[View online »](#)

**MW:** 42.6 kDa

**Gene Summary:** The protein encoded by this gene contains a RING zinc finger, a motif known to be involved in protein-protein interactions. The specific function of this gene has not yet been determined. Alternatively spliced transcript variants that encode the same protein have been reported. A pseudogene, which is also located on chromosome 3, has been defined for this gene. [provided by RefSeq, Jul 2008]