

## Product datasheet for **RC217197L3V**

### **PDE8A (NM\_173454) Human Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	PDE8A (NM_173454) Human Tagged ORF Clone Lentiviral Particle
Symbol:	PDE8A
Synonyms:	HsT19550
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_173454
ORF Size:	2349 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC217197).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_173454.1</a> , <a href="#">NP_775656.1</a>
RefSeq Size:	3793 bp
RefSeq ORF:	2352 bp
Locus ID:	5151
UniProt ID:	<a href="#">O60658</a>
Cytogenetics:	15q25.3
Protein Families:	Druggable Genome
Protein Pathways:	Progesterone-mediated oocyte maturation, Purine metabolism



[View online »](#)

**MW:** 88.1 kDa

**Gene Summary:** The protein encoded by this gene belongs to the cyclic nucleotide phosphodiesterase (PDE) family, and PDE8 subfamily. This PDE hydrolyzes the second messenger, cAMP, which is a regulator and mediator of a number of cellular responses to extracellular signals. Thus, by regulating the cellular concentration of cAMP, this protein plays a key role in many important physiological processes. Alternatively spliced transcript variants encoding different isoforms have been found for this gene.[provided by RefSeq, Jul 2011]