

Product datasheet for **RC217199L3V**

PDE9A (NM_001001581) Human Tagged ORF Clone Lentiviral Particle

Product data:

| | |
|---------------------------|--|
| Product Type: | Lentiviral Particles |
| Product Name: | PDE9A (NM_001001581) Human Tagged ORF Clone Lentiviral Particle |
| Symbol: | PDE9A |
| Synonyms: | HSPDE9A2 |
| Mammalian Cell Selection: | Puromycin |
| Vector: | pLenti-C-Myc-DDK-P2A-Puro (PS100092) |
| Tag: | Myc-DDK |
| ACCN: | NM_001001581 |
| ORF Size: | 1578 bp |
| ORF Nucleotide Sequence: | The ORF insert of this clone is exactly the same as(RC217199). |
| OTI Disclaimer: | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info |
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| RefSeq: | NM_001001581.1 , NP_001001581.1 |
| RefSeq Size: | 1954 bp |
| RefSeq ORF: | 1581 bp |
| Locus ID: | 5152 |
| UniProt ID: | O76083 |
| Cytogenetics: | 21q22.3 |
| Protein Families: | Druggable Genome |
| Protein Pathways: | Progesterone-mediated oocyte maturation, Purine metabolism |



[View online »](#)

MW: 61 kDa

Gene Summary: The protein encoded by this gene catalyzes the hydrolysis of cAMP and cGMP to their corresponding monophosphates. The encoded protein plays a role in signal transduction by regulating the intracellular concentration of these cyclic nucleotides. Multiple transcript variants encoding several different isoforms have been found for this gene. [provided by RefSeq, Jul 2008]