

## Product datasheet for **RC217202L4V**

### SLC2A11 (NM\_001024939) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	SLC2A11 (NM_001024939) Human Tagged ORF Clone Lentiviral Particle
Symbol:	SLC2A11
Synonyms:	GLUT10; GLUT11
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001024939
ORF Size:	1497 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC217202).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001024939.1</a>
RefSeq Size:	1590 bp
RefSeq ORF:	1500 bp
Locus ID:	66035
UniProt ID:	<a href="#">Q9BYW1</a>
Cytogenetics:	22q11.23
Protein Families:	Transmembrane
MW:	53.8 kDa



[View online »](#)

**Gene Summary:**

This gene belongs to a family of proteins that mediate the transport of sugars across the cell membrane. The encoded protein transports glucose and fructose. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Sep 2013]