

Product datasheet for **RC217216L4V**

Angiopoietin 2 (ANGPT2) (NM_001147) Human Tagged ORF Clone Lentiviral Particle

Product data:

| | |
|---------------------------|--|
| Product Type: | Lentiviral Particles |
| Product Name: | Angiopoietin 2 (ANGPT2) (NM_001147) Human Tagged ORF Clone Lentiviral Particle |
| Symbol: | ANGPT2 |
| Synonyms: | AGPT2; ANG2; LMPHM10 |
| Mammalian Cell Selection: | Puromycin |
| Vector: | pLenti-C-mGFP-P2A-Puro (PS100093) |
| Tag: | mGFP |
| ACCN: | NM_001147 |
| ORF Size: | 1488 bp |
| ORF Nucleotide Sequence: | The ORF insert of this clone is exactly the same as(RC217216). |
| OTI Disclaimer: | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info |
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| RefSeq: | NM_001147.1 |
| RefSeq Size: | 2269 bp |
| RefSeq ORF: | 1491 bp |
| Locus ID: | 285 |
| UniProt ID: | O15123 |
| Cytogenetics: | 8p23.1 |
| Domains: | FBG |
| Protein Families: | Druggable Genome, Secreted Protein |



[View online »](#)

MW: 56.95 kDa

Gene Summary: This gene belongs to the angiopoietin family of growth factors. The protein encoded by this gene is an antagonist of angiopoietin 1, and both angiopoietin 1 and angiopoietin 2 are ligands for the endothelial TEK receptor tyrosine kinase. Angiopoietin 2 is upregulated in multiple inflammatory diseases and is implicated in the direct control of inflammation-related signaling pathways. The encoded protein affects angiogenesis during embryogenesis and tumorigenesis, disrupts the vascular remodeling ability of angiopoietin 1, and may induce endothelial cell apoptosis. This gene serves a prognostic biomarker for acute respiratory distress syndrome. [provided by RefSeq, Aug 2020]