

Product datasheet for **RC217218L3V**

RCOR3 (NM_018254) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	RCOR3 (NM_018254) Human Tagged ORF Clone Lentiviral Particle
Symbol:	RCOR3
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_018254
ORF Size:	1485 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC217218).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_018254.2
RefSeq Size:	4083 bp
RefSeq ORF:	1488 bp
Locus ID:	55758
UniProt ID:	Q9P2K3
Cytogenetics:	1q32.2-q32.3
Domains:	ELM2, myb_DNA-binding
Protein Families:	Transcription Factors
MW:	55.4 kDa



[View online »](#)

Gene Summary:

May act as a component of a corepressor complex that represses transcription.
[UniProtKB/Swiss-Prot Function]