

Product datasheet for **RC217220L4V**

KHSRP (NM_003685) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	KHSRP (NM_003685) Human Tagged ORF Clone Lentiviral Particle
Symbol:	KHSRP
Synonyms:	FBP2; FUBP2; KSRP; p75
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_003685
ORF Size:	2133 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC217220).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_003685.1
RefSeq Size:	3009 bp
RefSeq ORF:	2136 bp
Locus ID:	8570
UniProt ID:	Q92945
Cytogenetics:	19p13.3
Protein Families:	Druggable Genome, Transcription Factors
MW:	72.9 kDa



[View online »](#)

Gene Summary:

The KHSRP gene encodes a multifunctional RNA-binding protein implicated in a variety of cellular processes, including transcription, alternative pre-mRNA splicing, and mRNA localization (Min et al., 1997 [PubMed 9136930]; Gherzi et al., 2004 [PubMed 15175153]). [supplied by OMIM, Apr 2010]