

Product datasheet for **RC217222A1V**

Human TMEM190 (NM_139172) AAV Particle

Product data:

Product Type: AAV Particles

Product Name: Human TMEM190 (NM_139172) AAV Particle

Tag: Myc-DDK

Symbol: TMEM190

Synonyms: MDAC1

Mammalian Cell Selection: None

Vector: pAAV-AC-Myc-DDK (PS100089)

ORF Nucleotide Sequence: >RC217222 ORF sequence
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGTTGGGCTGTGGGATCCCAGCGCTGGGCCTGCTCCTGCTGCTGCAGGGCTCGGCAGACGGAAATGGAA
TCCAGGGATTCTTCTACCCATGGAGCTGTGAGGGTGACATATGGGACCGGGAGAGCTGTGGGGGCCAGGC
GGCCATCGATAGCCCCAACCTCTGCCTGCGTCTCCGGTGTCTACCGCAATGGGGTCTGCTACCAACAG
CGTCCAGACGAAAACGTGCGGAGGAAGCACATGTGGGCGCTGGTCTGGACGTGCAGCGGCCTCCTCCTCC
TGAGTGCAGCATCTGCTTGTCTGGTGGGCCAAGCGCCGGGACGTGCTGCATATGCCCGTTTCTCTGGC
GGGTCCGTGTGACATGTCCAAGTCCGTCTCGTGTCTCCAAGCACCGAGGGACCAAGAAGACGCCGTCC
ACGGGCAGCGTGCCAGTCGCCCTGTCCAAGAGTCCAGGGATGTGGAGGAGGCACCGAGGGGAAGGGA
CGGAGGAGGGTGAGGAGACAGAGGGCGAGGAAGAGGAGGAT

ACGCGTACGCGGCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC217222 protein sequence
Red=Cloning site Green=Tags(s)

MLGCGIPALGLLLLQGSADGNGIQGFFYPWSEGDIDWRESCGGQAAIDSPNLCLRLRCCYRNGVCYHQ
RPDENVRKHMWALVWTCGLLLLSCSICLFWWAKRRDVLHMPGLAGPCDMSKSVSLLSKHRGTTKTPS
TGSVPVALSKESRDVEGGTEGEGTEEGEETEGEEED

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species: Human

Serotype: AAV-2



[View online »](#)

ACCN:	NM_139172
ORF Size:	531 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_139172.1</u>
RefSeq Size:	593 bp
RefSeq ORF:	534 bp
Locus ID:	147744
UniProt ID:	<u>Q8WZ59</u>
Cytogenetics:	19q13.42
MW:	19.5 kDa