

Product datasheet for **RC217226L4V**

HOXC5 (NM_018953) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	HOXC5 (NM_018953) Human Tagged ORF Clone Lentiviral Particle
Symbol:	HOXC5
Synonyms:	CP11; HOX3; HOX3D
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_018953
ORF Size:	666 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC217226).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_018953.2 , NP_061826.1
RefSeq Size:	1613 bp
RefSeq ORF:	669 bp
Locus ID:	3222
UniProt ID:	Q00444
Cytogenetics:	12q13.13
Domains:	homeobox
Protein Families:	Transcription Factors



[View online »](#)

MW: 24.8 kDa

Gene Summary: This gene belongs to the homeobox family of genes. The homeobox genes encode a highly conserved family of transcription factors that play an important role in morphogenesis in all multicellular organisms. Mammals possess four similar homeobox gene clusters, HOXA, HOXB, HOXC and HOXD, which are located on different chromosomes and consist of 9 to 11 genes arranged in tandem. This gene, HOXC5, is one of several homeobox HOXC genes located in a cluster on chromosome 12. Three genes, HOXC5, HOXC4 and HOXC6, share a 5' non-coding exon. Transcripts may include the shared exon spliced to the gene-specific exons, or they may include only the gene-specific exons. Two alternatively spliced variants have been described for HOXC5. The transcript variant which includes the shared exon apparently doesn't encode a protein. The protein-coding transcript variant contains gene-specific exons only. [provided by RefSeq, Jul 2008]