

Product datasheet for **RC217232**

Choline Acetyltransferase (CHAT) (NM_020986) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Choline Acetyltransferase (CHAT) (NM_020986) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	CHAT
Synonyms:	CHOACTASE; CMS1A; CMS1A2; CMS6
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin



[View online »](#)

ORF Nucleotide Sequence:

>RC217232 representing NM_020986
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**GCGATCGCC**

ATGGCAGCAAAAACCTCCAGCAGTGTAGGAGTCTGGGCTGCCCAAACCTGCCCGTGCCCGCTGCAGCAGA
 CCCTGGCCACGTACCTGCAGTGCATGCGACACTTGGTGTCTGAGGAGCAGTTTCAGGAAGAGCCAGGCCAT
 TGTGCAGCAGTTTGGGGCCCTGGTGGCCTCGGCGAGACCCTGCAGCAGAAAACCTCGGAGCGGCAGGAG
 AAGACAGCCAACCTGGGTGTCTGAGTACTGGCTGAATGACATGTATCTCAACAACCGCCTGGCCCTGCCTG
 TCAACTCCAGCCCTGCCGTGATCTTTGCTCGGCAGCACTTCCCTGGCACCAGTACCAGCTGAGGTTTGC
 AGCCAGCCTCATCTCTGGTGTACTCAGCTACAAGGCCCTGCTGGACAGCCACTCCATTCCCAGTACTGT
 GCCAAAGGCCAGCTGTCAGGGCAGCCCTTTGCATGAAGCAATACTATGGGCTCTTCTCCTCTACCGGC
 TCCCGGCCATACCCAGGACACGCTGGTGGCTCAGAACAGCAGCATCATGCCGGAGCCTGAGCACGTCT
 CGTAGCCTGCTCAATCAGTTCTTTGTCTTGGATGTTGTCATTAATTTCCGCCGTCTCAGTGTAGGGGGAT
 CTGTCTACTCAGTTGAGAAAAGATAGTCAAAATGGCTTCCAACGAGGACGAGCGTTTGCCTCCAATTGGCC
 TGCTGACGTCTGACGGGAGGAGCGAGTGGGCCGAGGCCAGGACGGTCTCGTGAAAGACTCCACCAACCG
 GGACTCGCTGGACATGATTGAGCGCTGCATCTGCCTTGTATGCCTGGACGCGCCAGGAGCGTGGAGCTC
 AGCGACACCCACAGGGCACTCCAGCTCCTTACGGCGGAGGCTACAGCAAGAACGGGCAATCGCTGGT
 ACGACAAGTCCCTGCAGTTTGTGGTGGGCCGAGACGGCACCTGCGGTGTGGTGTGCGAACACTCCCCATT
 CGATGGCATCGTCTGGTGCAGTGCAGTGCATCTGCTCAAGCACATGACGAGAGCAGCAGGAAGCTG
 ATCCGAGCAGACTCCGTCAGCGAGCTCCCCGCCCGGGAGGCTGCGGTGAAATGCTCCCCGAAATTC
 AAGGCCACTTAGCCTCCTCGGCAGAAAACCTTCAACGAATAGTAAAGAACCTTGACTTCATTGTCTATAA
 GTTTGACAACTATGGGAAAACATTCATTAAGAAGCAGAAAATGCAGCCCTGATGCCTTCATCCAGGTGCC
 CTCCAGCTGGCCTTCTACAGGCTCCATCGAAGACTGGTGCCACCTACGAGAGCGCTCCATCCGCCGAT
 TCCAGGAGGGACGCGTGGACAACATCAGATCGGCCACTCCAGAGGCACTGGCTTTTGTGAGAGCCGTGAC
 TGACCACAAGGCTGCTGTGCCAGCTTCTGAGAAGCTTCTGCTCCTGAAGGATGCCATCCGTGCCAGACT
 GCATACACAGTGCATGGCCATAACAGGGATGGCCATTGACAACCACCTGCTGGCACTGCGGGAGCTGGCC
 GGGCCATGTGCAAGGAGCTGCCCGAGATGTTTCATGGATGAAACCTACCTGATGAGCAACCGGTTTGTCT
 CTCCTACTGCCAGGTGCCACAACCACGGAGATGTTCTGCTGCTATGGTCTGTGGTCCCAAATGGGTAT
 GGTGCCCTGCTACAACCCCGCCAGAGACCATCTTTTCTGCATCTCTAGCTTTCACAGCTGCAAAGAGA
 CTCTTTCTAGCAAGTTTGCAAAAGCTGTGGAAGAAAGCCTCATTGACATGAGAGACCTCTGCAGTCTGCT
 GCCGCTACTGAGAGCAAGCCATTGGCAACAAAGGAAAAGCCACGAGGCCAGCCAGGGACACCAACCT

AG**GCGACCG**ACGCGTACGCGGCCGCTCGAGCAGAAAACCTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCC
 TGGATTACAAGGATGACGACGATAAGGTTTAA

Protein Sequence:


>RC217232 representing NM_020986
 Red=Cloning site Green=Tags(s)

MAAKTPSSEESGLPKLPVPLQQTALATYLQCMRHLVSEEQFRKSQAIIVQFQFAGPGLGETLQQLLERQE
 KTANWVSEYWLNDMYLNNRLALPVNSSPAVIFARQHFPGTDDQLRFAASLISGLVSYKALLDSHSIPTDC
 AKGQLSGQPLCMKQYYGLFSSYRLPGHTQDTLVAQNSSIMPEPEHVIVACCNQFFVLDVVINFRRLESEG
 LFTQLRKIVKMASNEDERLPIIGLLTSDGRSEWAEARTVLVKDSTNRDSLDMIERICLVCLDAPGGVEL
 SDTHRALQLLHGGYSKNGANRWYDKSLQFVVGRTDGTGCVVCEHSPFDGIVLVQTEHLLKHMTQSSRKL
 IRADSVSELPAAPRRLRWKCSPEIQGHLASSAEKLRIVKNLDFIVYKFDNYGKTFIKKQKSPDAFIQVA
 LQLAFYRLHRRLVPTYESASIRRFQEGVRDNIIRSATPEALAFVRAVTDHKAAPASEKLLLLKDAIRAQT
 AYTVMAITGMAIDNHLALRELARAMCKELPEMFMDETYLSNRFLVSTSQVPTTTEMFCYGPVVPNGY
 GACYNPQPETILFCISSFHCKETSSSKFAKAVEESLIDMRDLCSLLPPTESKPLATKEKATRPSQGHQP


SGPTRTRRLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-RsrII

Cloning Scheme:



Plasmid Map:



ACCN: NM_020986

ORF Size: 1890 bp

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_020986.3 , NP_066266.3
RefSeq Size:	2114 bp
RefSeq ORF:	1893 bp
Locus ID:	1103
UniProt ID:	P28329
Cytogenetics:	10q11.23
Protein Families:	Druggable Genome
Protein Pathways:	Glycerophospholipid metabolism
MW:	70.4 kDa
Gene Summary:	This gene encodes an enzyme which catalyzes the biosynthesis of the neurotransmitter acetylcholine. This gene product is a characteristic feature of cholinergic neurons, and changes in these neurons may explain some of the symptoms of Alzheimer's disease. Polymorphisms in this gene have been associated with Alzheimer's disease and mild cognitive impairment. Mutations in this gene are associated with congenital myasthenic syndrome associated with episodic apnea. Multiple transcript variants encoding different isoforms have been found for this gene, and some of these variants have been shown to encode more than one isoform. [provided by RefSeq, May 2010]