

Product datasheet for RC217241

SPLUNC3 (BPIFA3) (NM_178466) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	SPLUNC3 (BPIFA3) (NM_178466) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	SPLUNC3
Synonyms:	C20orf71; SPLUNC3
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>RC217241 representing NM_178466 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGATGTGTCCACTCTGGAGGCTCCTCATCTTCTCGGGTTGCTGGCCTTGCCCTTGGCACCACACAAGC
AGCCTTGGCCTGGCCTGGCCCAAGCCCACAGAGACAACAAATCCACCCTGGCAAGAATTATTGCTCAGGG
CCTCATAAAGCACAAACGCAGAAAGCCGAATTCAGAACATCCACTTTGGGGACAGACTGAATGCCTCAGCA
CAAGTGGCCCCAGGGCTGGTGGCTGGTAATCAGCGGCAGGAAACACCAGCAGCAAGAGAGCAGCA
TCAACATCACCAACATTCAGCTGGACTGTGGTGGGATCCAGATATCATTCCATAAGGAGTGGTTCTCGGC
AAATATCTCACTTGAATTTGACCTTGAATTGAGACCGTCCTTCGATAACAACATCGTAAAGATGTGTGCA
CATATGAGCATCGTTGTGGAGTTCTGGCTGGAGAAAGACGAGTTTGGCCGGAGGGATCTGGTGTATAGGCA
AATGCGATGCAGAGCCCAGCAGTGTCCATGTGGCCATCCTCACTGAGGCTATCCCACAAAGATGAATCA
GTTTCTCTACAACCTCAAAGAGAATCTGCAAAAAGTTCTCCCACACATGGTAGAAAGTCAGGTATGTCTC
CTGATCGGTGAAATCCTCGGCAGCTGGATGTGAAACTGTTGAAAAGCCTCATAGAACAGGAGGCTGCTC
ATGAACCAACCCACCATGAAACCAGCCAACCCTCTGCATGCCAGGCTGGAGAGTCCCCCAGC

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_178466.5](#)

RefSeq Size: 1112 bp

RefSeq ORF: 765 bp


Locus ID: 128861

UniProt ID: [Q9BQP9](#)

Cytogenetics: 20q11.21

MW: 28.3 kDa

Product images:



Circular map for RC217241