

Product datasheet for **RC217289A1V**

Human TIM 1 (HAVCR1) (NM_001099414) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human TIM 1 (HAVCR1) (NM_001099414) AAV Particle
Tag:	Myc-DDK
Symbol:	TIM 1
Synonyms:	HAVCR; HAVCR-1; KIM-1; KIM1; TIM; TIM-1; TIM1; TIMD-1; TIMD1
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC217289 representing NM_001099414 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGC**

ATGCATCCTCAAGTGGTCATCTTAAGCCTCATCTACATCTGGCAGATTCTGTAGCTGGTTCTGTAAAGG
TTGGTGGAGAGGCAGGTCCATCTGTACACTACCCTGCCACTACAGTGGAGCTGTACATCCATGTGCTG
GAATAGAGGCTCATGTTCTCTATTACATGCCAAATGGCATTGTCTGGACCAATGGAACCCACGTCACC
TATCGGAAGGACACACGCTATAAGCTATTGGGGGACCTTTCAAGAAGGGATGTCTCTTTGACCATAGAAA
ATACAGCTGTGTCTGACAGTGGCGTATATTGTTGCCGTGTTGAGCACCGTGGTGGTTCAATGACATGAA
AATCACCGTATCATTGGAGATTGTGCCACCAAGGTCACGACTACTCCAATTGTCACAACTGTTCCAACC
GTCACGACTGTTCAACGAGCACCCTGTTCCAACGACAACGACTGTTCCAATGACGACTGTTCCAACGA
CAACTGTTCCAACAACATGAGCATTCCAACGACAACGACTGTTCTGACGACAATGACTGTTCAACGAC
AACGAGCGTTCCAACGACAACGAGCATTCCAACAACAAGTGTCCAGTGACAACAACGACTGTTCTTACC
TTTGTTCCTCCAATGCCTTTGCCAGGCAGAACCATGAACAGTAGCCACTTCACCATCTTCACCTCAGC
CAGCAGAAACCCACCCTACGACTGCAGGGAGCAATAAGGAGAGAACCACACGCTCACCATTGTAATC
TTACACAACAGATGGGAATGACACCGTGACAGAGTCTTCAGATGGCCTTTGGAATAACAATCAAACCTCA
CTGTTCTAGAACATAGTCTACTGACGGCAATACCACTAAAGGAATCTATGCTGGAGTCTGTATTCTG
TCTTGGTGCTTCTTGCTCTTTGGGTGTCTATTGCCAAAAAGTATTTCTTCAAAAAGGAGGTTCAACA
ACTAAGTGTTCATTTAGCAGCCTTCAAATTAAGCTTTGCAAAATGCAGTTGAAAAGGAAGTCCAAGCA
GAAGACAATATCTACATTGAGAATAGTCTTTATGCCACGGAC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC217289 representing NM_001099414 Red=Cloning site Green=Tags(s) MHPQVVILSLILHLADSVAGSVKVGGEAGPSVTLPCHYSGAVTSMCWNRGSCSLFTCQNGI VWTNGTHVT YRKDTRYKLLGDL SRRDVSLTIENTAVSDSGVYCCRVEHRGW FNDMKITVSLEIVPPKVTTTPIVTTVPT VTTVRTSTTVPTTTVPMTTVPTTVPITMSIPTTTTTLTMTVSTTTSVPTTTSIPTTTSVPVTTTVST FVPPMPLPRQNHVPATSPSSPQPAETHPTTLQGAIRREPTSSPLYSYTTDGNDTVTESSDGLWNNQTQ LFLEHSLLTANTTKGIYAGVCISVLVLLALLGVIIAKKYFFKKEVQQLSVSFSSSLQIKALQNAVEKEVQA EDNIYIENSLYATD TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_001099414
ORF Size:	1092 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_001099414.1</u> , <u>NP_001092884.1</u>
RefSeq Size:	1493 bp
RefSeq ORF:	1095 bp
Locus ID:	26762
Cytogenetics:	5q33.3
MW:	39.25 kDa