

Product datasheet for RC217291L1

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

OriGene Technologies, Inc.

GLI2 (NM_005270) Human Tagged Lenti ORF Clone

Product data:

Product Type: Expression Plasmids

Product Name: GLI2 (NM_005270) Human Tagged Lenti ORF Clone

Tag: Myc-DDK

Symbol: GLI2

Synonyms: CJS; HPE9; PHS2; THP1; THP2

Mammalian Cell None

Selection:

Vector:pLenti-C-Myc-DDK (PS100064)E. coli Selection:Chloramphenicol (34 ug/mL)

ORF Nucleotide T


JRF Nucleotide


The ORF insert of this clone is exactly the same as(RC217291).

Sequence:

Restriction Sites: Sgfl-Notl

Cloning Scheme:





^{*} The last codon before the Stop codon of the ORF.

ACCN: NM_005270

ORF Size: 4758 bp



GLI2 (NM_005270) Human Tagged Lenti ORF Clone - RC217291L1

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of

reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing

variants is recommended prior to use. More info

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression

varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube

containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method: 1. Centrifuge at 5,000xg for 5min.

2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.

3. Close the tube and incubate for 10 minutes at room temperature.

4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid

at the bottom.

5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of

shipping when stored at -20°C.

RefSeq: <u>NM 005270.3</u>

RefSeq Size: 6549 bp
RefSeq ORF: 4761 bp
Locus ID: 2736

UniProt ID: P10070
Cytogenetics: 2q14.2

Protein Families: Adult stem cells, Cancer stem cells, Druggable Genome, Embryonic stem cells, ES Cell

Differentiation/IPS

Protein Pathways: Basal cell carcinoma, Hedgehog signaling pathway, Pathways in cancer

MW: 167.6 kDa

Gene Summary: This gene encodes a protein which belongs to the C2H2-type zinc finger protein subclass of

the Gli family. Members of this subclass are characterized as transcription factors which bind DNA through zinc finger motifs. These motifs contain conserved H-C links. Gli family zinc finger proteins are mediators of Sonic hedgehog (Shh) signaling and they are implicated as potent oncogenes in the embryonal carcinoma cell. The protein encoded by this gene localizes to the cytoplasm and activates patched Drosophila homolog (PTCH) gene expression. It is also thought to play a role during embryogenesis. The encoded protein is

associated with several phenotypes- Greig cephalopolysyndactyly syndrome, Pallister-Hall syndrome, preaxial polydactyly type IV, postaxial polydactyly types A1 and B. [provided by

RefSeq, Jul 2008]



Product images:

