

## Product datasheet for **RC217297L1V**

### **CPT1B (NM\_004377) Human Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	CPT1B (NM_004377) Human Tagged ORF Clone Lentiviral Particle
Symbol:	CPT1B
Synonyms:	CPT1-M; CPT1M; CPTI; CPTI-M; M-CPT1; MCCPT1; MCPT1
Mammalian Cell Selection:	None
Vector:	pLenti-C-Myc-DDK (PS100064)
Tag:	Myc-DDK
ACCN:	NM_004377
ORF Size:	2316 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC217297).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_004377.2</a>
RefSeq Size:	2879 bp
RefSeq ORF:	2319 bp
Locus ID:	1375
UniProt ID:	<a href="#">Q92523</a>
Cytogenetics:	22q13.33
Domains:	Carn_acyltransf
Protein Families:	Druggable Genome, Transmembrane



[View online »](#)

**Protein Pathways:** Adipocytokine signaling pathway, Fatty acid metabolism, PPAR signaling pathway

**MW:** 87.6 kDa

**Gene Summary:** The protein encoded by this gene, a member of the carnitine/choline acetyltransferase family, is the rate-controlling enzyme of the long-chain fatty acid beta-oxidation pathway in muscle mitochondria. This enzyme is required for the net transport of long-chain fatty acyl-CoAs from the cytoplasm into the mitochondria. Multiple transcript variants encoding different isoforms have been found for this gene, and read-through transcripts are expressed from the upstream locus that include exons from this gene. [provided by RefSeq, Jun 2009]