

Product datasheet for **RC217354L2V**

DNA polymerase eta (POLH) (NM_006502) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	DNA polymerase eta (POLH) (NM_006502) Human Tagged ORF Clone Lentiviral Particle
Symbol:	DNA polymerase eta
Synonyms:	RAD30; RAD30A; XP-V; XPV
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_006502
ORF Size:	2139 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC217354).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_006502.1
RefSeq Size:	3464 bp
RefSeq ORF:	2142 bp
Locus ID:	5429
UniProt ID:	Q9Y253
Cytogenetics:	6p21.1
Domains:	IMS
Protein Families:	Druggable Genome


[View online »](#)

MW: 78.2 kDa

Gene Summary: This gene encodes a member of the Y family of specialized DNA polymerases. It copies undamaged DNA with a lower fidelity than other DNA-directed polymerases. However, it accurately replicates UV-damaged DNA; when thymine dimers are present, this polymerase inserts the complementary nucleotides in the newly synthesized DNA, thereby bypassing the lesion and suppressing the mutagenic effect of UV-induced DNA damage. This polymerase is thought to be involved in hypermutation during immunoglobulin class switch recombination. Mutations in this gene result in XPV, a variant type of xeroderma pigmentosum. Several transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, May 2014]