

## Product datasheet for **RC217360L4V**

### CDw75 (ST6GAL1) (NM\_173217) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	CDw75 (ST6GAL1) (NM_173217) Human Tagged ORF Clone Lentiviral Particle
Symbol:	CDw75
Synonyms:	SIAT1; ST6Gall; ST6N
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_173217
ORF Size:	525 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC217360).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_173217.1</a>
RefSeq Size:	3746 bp
RefSeq ORF:	528 bp
Locus ID:	6480
UniProt ID:	<a href="#">P15907</a>
Cytogenetics:	3q27.3
Protein Families:	Secreted Protein
Protein Pathways:	Metabolic pathways, N-Glycan biosynthesis



[View online »](#)

**MW:** 20.6 kDa

**Gene Summary:** This gene encodes a member of glycosyltransferase family 29. The encoded protein is a type II membrane protein that catalyzes the transfer of sialic acid from CMP-sialic acid to galactose-containing substrates. The protein, which is normally found in the Golgi but can be proteolytically processed to a soluble form, is involved in the generation of the cell-surface carbohydrate determinants and differentiation antigens HB-6, CD75, and CD76. This gene has been incorrectly referred to as CD75. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jul 2017]