

Product datasheet for **RC217371L4V**

PTPN5 (NM_032781) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	PTPN5 (NM_032781) Human Tagged ORF Clone Lentiviral Particle
Symbol:	PTPN5
Synonyms:	PTPSTEP; STEP; STEP61
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_032781
ORF Size:	1695 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC217371).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_032781.2
RefSeq Size:	2965 bp
RefSeq ORF:	1698 bp
Locus ID:	84867
UniProt ID:	P54829
Cytogenetics:	11p15.1
Domains:	Y_phosphatase, PTPc_motif
Protein Families:	Druggable Genome, Phosphatase, Transmembrane



[View online »](#)

Protein Pathways: MAPK signaling pathway

MW: 63.4 kDa

Gene Summary: May regulate the activity of several effector molecules involved in synaptic plasticity and neuronal cell survival, including MAPKs, Src family kinases and NMDA receptors.
[UniProtKB/Swiss-Prot Function]