

## Product datasheet for **RC217381L2V**

### **FAM36A (COX20) (NM\_198076) Human Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	FAM36A (COX20) (NM_198076) Human Tagged ORF Clone Lentiviral Particle
Symbol:	FAM36A
Synonyms:	FAM36A; MC4DN11
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_198076
ORF Size:	354 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC217381).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_198076.4</a>
RefSeq Size:	2632 bp
RefSeq ORF:	357 bp
Locus ID:	116228
UniProt ID:	<a href="#">Q5RI15</a>
Cytogenetics:	1q44
MW:	13.3 kDa



[View online »](#)

**Gene Summary:**

This gene encodes a protein that plays a role in the assembly of cytochrome C oxidase, an important component of the respiratory pathway. It contains two transmembrane helices and localizes to the mitochondrial membrane. Mutations in this gene can cause mitochondrial complex IV deficiency, which results in ataxia and muscle hypotonia. There are multiple pseudogenes for this gene. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Aug 2015]