

Product datasheet for **RC217405L4V**

USP42 (NM_032172) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	USP42 (NM_032172) Human Tagged ORF Clone Lentiviral Particle
Symbol:	USP42
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_032172
ORF Size:	3948 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC217405).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_032172.1
RefSeq Size:	5109 bp
RefSeq ORF:	3951 bp
Locus ID:	84132
UniProt ID:	Q9H9J4
Cytogenetics:	7p22.1
Protein Families:	Druggable Genome, Protease
MW:	144.2 kDa



[View online »](#)

Gene Summary: Deubiquitinating enzyme which may play an important role during spermatogenesis.
[UniProtKB/Swiss-Prot Function]