

Product datasheet for **RC217406L1V**

VGLL2 (NM_182645) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	VGLL2 (NM_182645) Human Tagged ORF Clone Lentiviral Particle
Symbol:	VGLL2
Synonyms:	VGL2; VITO1
Mammalian Cell Selection:	None
Vector:	pLenti-C-Myc-DDK (PS100064)
Tag:	Myc-DDK
ACCN:	NM_182645
ORF Size:	951 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC217406).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_182645.2 , NP_872586.1
RefSeq Size:	2212 bp
RefSeq ORF:	954 bp
Locus ID:	245806
UniProt ID:	Q8N8G2
Cytogenetics:	6q22.1
MW:	33.2 kDa



[View online »](#)

Gene Summary:

This gene encodes a protein with a transcriptional enhancer factor 1 (TEF-1) interaction domain. The encoded protein may act as a co-factor of TEF-1 regulated gene expression during skeletal muscle development. Alternatively spliced transcript variants encoding multiple isoforms have been observed. [provided by RefSeq, Jul 2012]