

Product datasheet for **RC217412L1V**

PMS1 (NM_000534) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	PMS1 (NM_000534) Human Tagged ORF Clone Lentiviral Particle
Symbol:	PMS1
Synonyms:	HNPCC3; hPMS1; MLH2; PMSL1
Mammalian Cell Selection:	None
Vector:	pLenti-C-Myc-DDK (PS100064)
Tag:	Myc-DDK
ACCN:	NM_000534
ORF Size:	2796 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC217412).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_000534.3
RefSeq Size:	3239 bp
RefSeq ORF:	2799 bp
Locus ID:	5378
UniProt ID:	P54277
Cytogenetics:	2q32.2
Domains:	HMG, DNA_mis_repair, HATPase_c
Protein Families:	Druggable Genome, Transcription Factors



[View online »](#)

MW: 105.7 kDa

Gene Summary: This gene encodes a protein belonging to the DNA mismatch repair mutL/hexB family. This protein is thought to be involved in the repair of DNA mismatches, and it can form heterodimers with MLH1, a known DNA mismatch repair protein. Mutations in this gene cause hereditary nonpolyposis colorectal cancer type 3 (HNPCC3) either alone or in combination with mutations in other genes involved in the HNPCC phenotype, which is also known as Lynch syndrome. [provided by RefSeq, Jul 2008]