

Product datasheet for **RC217447L1V**

HFE (NM_000410) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	HFE (NM_000410) Human Tagged ORF Clone Lentiviral Particle
Symbol:	HFE
Synonyms:	HFE1; HH; HLA-H; MVCD7; TFQTL2
Mammalian Cell Selection:	None
Vector:	pLenti-C-Myc-DDK (PS100064)
Tag:	Myc-DDK
ACCN:	NM_000410
ORF Size:	1044 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC217447).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_000410.2
RefSeq Size:	2727 bp
RefSeq ORF:	1047 bp
Locus ID:	3077
UniProt ID:	Q30201
Cytogenetics:	6p22.2
Domains:	MHC_I, ig, IGc1
Protein Families:	Druggable Genome, Transmembrane




[View online »](#)

MW: 40.11 kDa

Gene Summary: The protein encoded by this gene is a membrane protein that is similar to MHC class I-type proteins and associates with beta2-microglobulin (beta2M). It is thought that this protein functions to regulate iron absorption by regulating the interaction of the transferrin receptor with transferrin. The iron storage disorder, hereditary haemochromatosis, is a recessive genetic disorder that results from defects in this gene. At least nine alternatively spliced variants have been described for this gene. Additional variants have been found but their full-length nature has not been determined. [provided by RefSeq, Jul 2008]

Product images:



[RC217447L1] was used to prepare Lentiviral particles using [TR30037] packaging kit. HEK293T cells were transduced with RC217447L1V particle to overexpress human HFE-Myc-DDK fusion protein.