

## Product datasheet for **RC217455L3V**

### **CARD14 (NM\_024110) Human Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	CARD14 (NM_024110) Human Tagged ORF Clone Lentiviral Particle
Symbol:	CARD14
Synonyms:	BIMP2; CARMA2; PRP; PSORS2; PSS1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_024110
ORF Size:	3012 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC217455).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_024110.2</a>
RefSeq Size:	3917 bp
RefSeq ORF:	3015 bp
Locus ID:	79092
UniProt ID:	<a href="#">Q9BXL6</a>
Cytogenetics:	17q25.3
Protein Families:	Druggable Genome
MW:	113.1 kDa



[View online »](#)

**Gene Summary:**

This gene encodes a caspase recruitment domain-containing protein that is a member of the membrane-associated guanylate kinase (MAGUK) family of proteins. Members of this protein family are scaffold proteins that are involved in a diverse array of cellular processes including cellular adhesion, signal transduction and cell polarity control. This protein has been shown to specifically interact with BCL10, a protein known to function as a positive regulator of cell apoptosis and NF-kappaB activation. Alternate splicing results in multiple transcript variants. [provided by RefSeq, Apr 2012]