

Product datasheet for **RC217464L4V**

NETO1 (NM_138999) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	NETO1 (NM_138999) Human Tagged ORF Clone Lentiviral Particle
Symbol:	NETO1
Synonyms:	BCTL1; BTCL1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_138999
ORF Size:	468 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC217464).
OTI Disclaimer:	<p>Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.</p> <p>The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info</p>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_138999.1
RefSeq Size:	1849 bp
RefSeq ORF:	471 bp



[View online »](#)

Locus ID:	81832
UniProt ID:	Q8TDF5
Cytogenetics:	18q22.3
Protein Families:	Druggable Genome, Transmembrane
MW:	17.7 kDa
Gene Summary:	<p>This gene encodes a transmembrane protein containing two extracellular CUB domains followed by a low-density lipoprotein class A (LDLa) domain. This protein is thought to play a critical role in spatial learning and memory by regulating the function of synaptic N-methyl-D-aspartic acid receptor complexes in the hippocampus. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene. [provided by RefSeq, Aug 2017]</p>