

Product datasheet for **RC217503L3V**

ST2 (IL1RL1) (NM_016232) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	ST2 (IL1RL1) (NM_016232) Human Tagged ORF Clone Lentiviral Particle
Symbol:	IL1RL1
Synonyms:	DER4; FIT-1; IL33R; ST2; ST2L; ST2V; T1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_016232
ORF Size:	1668 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC217503).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_016232.4
RefSeq Size:	2058 bp
RefSeq ORF:	1671 bp
Locus ID:	9173
UniProt ID:	Q01638
Cytogenetics:	2q12.1
Domains:	TIR, ig, IGc2, IG
Protein Families:	Druggable Genome, Secreted Protein, Transmembrane




[View online »](#)

MW: 63.36 kDa

Gene Summary: The protein encoded by this gene is a member of the interleukin 1 receptor family. Studies of the similar gene in mouse suggested that this receptor can be induced by proinflammatory stimuli, and may be involved in the function of helper T cells. This gene, interleukin 1 receptor, type I (IL1R1), interleukin 1 receptor, type II (IL1R2) and interleukin 1 receptor-like 2 (IL1RL2) form a cytokine receptor gene cluster in a region mapped to chromosome 2q12. Alternative splicing of this gene results in multiple transcript variants. [provided by RefSeq, Jul 2008]

Product images:



[RC217503L3] was used to prepare Lentiviral particles using [TR30037] packaging kit. HEK293T cells were transduced with RC217503L3V particle to overexpress human IL1RL1-Myc-DDK fusion protein.